


The Use of NMR Spectroscopy for Chiral Discrimination

Thomas J. Wenzel
Department of Chemistry
Bates College
Lewiston, Maine



2007

“Using NMR Spectroscopic Methods to Determine Enantiomeric Purity and Assign Absolute Stereochemistry,” Wenzel, T.J.; Chisholm, C.D., *Progress in NMR Spectroscopy*, (DOI:10.1016/j.pnmrs.2010.07.003).

“Assignment of Absolute Configuration Using Chiral Reagents and NMR Spectroscopy,” Wenzel, T.J.; Chisholm, C.D., *Chirality*, (DOI: 10.1002/chir.20889).

Categories of Reagents

- Chiral Derivatizing Agents
- Chiral Solvating Agents
- Metal Complexes
- Liquid Crystals

Chiral Derivatizing Agents Raban and Mislow (1965)

- Form a covalent bond between an optically pure reagent ((*S*)-CDA) and the compound of interest (Sub)



- Resulting compounds are diastereomers
- Signals double in NMR spectrum (Chemical Shift Anisotropy) – areas proportional to percent of each enantiomer

Chiral Derivatizing Agents: Key Criteria if Using to Determine Enantiomeric Excess

- No racemization
- No kinetic resolution
- Need 100% enantiomeric purity of the reagent

Chiral Solvating Agents (Pirkle – 1966) Metal Complexes (Whitesides and Lewis – 1970)

- Form non-covalent interactions between an optically pure reagent ((*S*)-CSA) and the compound of interest (Sub)



- Resulting compounds are diastereomers
- K_R and K_S are likely different – causes different time-averaged solvation environments


Chiral Solvating Agents and Metal Complexes

- Mix directly in an NMR tube
- Preferable to have fast exchange – NMR spectrum is a time average of bound and unbound forms (CSA + Sub = CSA-Sub)
- High concentration of CSA usually leads to larger discrimination
- Often see enhanced enantiomeric discrimination at lower temperatures
- CSA does not need to be 100% enantiomerically pure


Assigning Absolute Stereochemistry


- Mechanism of discrimination is understood and characteristic changes in chemical shifts occur in the spectrum
 - More common with certain families of chiral derivatizing agents
 - Possible with some chiral solvating agents
- Empirical trend
 - Best if use known model compounds as close as possible in structural features to the unknown

Mosher Method: α -methoxy- α -trifluoromethylphenylacetic acid - MTPA (Dale and Mosher – 1973)



- Prepare derivatives with (*R*)- and (*S*)-forms of the reagent (esters of secondary alcohols)
- *Syn-periplanar* arrangement of HC-O-C(O)-C atoms (secondary alcohols)
- Calculate $\Delta\delta^{RS}$ values – negative for L_1 , positive for L_2






$\Delta\delta^{RS}$ depends on:

- Extent of conformational preference/how it influences the shielding
- Degree of shielding
(anthryl > naphthyl > phenyl)


2,2,2-Trifluoro-1-(9-anthryl)ethanol (TFAE) (Pirkle's Alcohol)



Versatile chiral solvating agent

- Can determine optical purity
- Can assign absolute configurations for certain classes of compounds

Absolute Configurations - TFAE



Sulfoxides

Metal Complexes: Expand Coordination Number or Displace Ligand (Donor/Acceptor Association)

- Lanthanides – Hard Lewis bases
 - Nitrogen- and Oxygen-containing compounds
- Platinum, Palladium, Rhodium and Silver – Soft Lewis bases
 - Alkenes, alkynes, aromatics, phosphorus-containing, sulfur-containing, alkyl halides

Liquid Crystals

Sackmann, Meiboom, Snyder (1968)


- Forms ordered material in a magnetic field
- Pair of enantiomers have different molecular orientations in the liquid crystal
- Three discrimination mechanisms
 - Chemical shift anisotropy (least useful)
 - Different dipolar coupling constants (^1H - ^{13}C)
 - Differences in quadrupolar splitting (^2H) (most useful)

Quadrupolar Splitting

- Not observed in solution because of rapid tumbling
- Observed in ordered media and extent of splitting depends on orientation relative to the applied magnetic field

2-²H-Propionic Acid

Proton-decoupled deuterium NMR spectrum



Poly(γ -benzyl-L-glutamate) – (PBLG)

Incredible Versatility


- Only need different packing orders
- Do not need specific interactions between the substrate and the liquid crystal
- Effective for virtually any class of compound
 - Includes aliphatic hydrocarbons
- Especially effective for resonances of nuclei remote to the chiral center

Deuterium Labeling

- Only need deuterium as a signal – better to use achiral reagents so no concern about kinetic resolution or racemization
 - Convert $-\text{CO}_2\text{H}$ to $-\text{CO}_2\text{CD}_3$
 - Add perdeutero benzoyl group (have o-, m- and p-protons as potential probes)
- Provides a single, strong signal (or a few easily assigned signals) for the analysis

Crown Ethers


(18-crown-6)-2,3,11,12-tetracarboxylic acid





Commercially Available



Association of Primary Amines




- Wenzel, T. J.; Freeman, B. E.; Sek, D. C.; Zopf, J. J.; Nakamura, T.; Yongzhu, J.; Hirose, K.; Tobe, Y, *Analytical and Bioanalytical Chemistry*, **2004**, *378*, 1536-1547.
- Wenzel, T. J.; Thurston, J. E., *Tetrahedron Letters*, **2000**, *41*, 3769-3772.
- Wenzel, T. J.; Thurston, J. E., *Journal of Organic Chemistry*, **2000**, *65*, 1243-1248.





Wenzel, T.J.; Bourne, C.E.; Clark, R.L., *Tetrahedron: Asymmetry*, **2009**, *20*, 2052-2060.
 Chisholm, C.D.; Fülöp, F.; Forró, E.; Wenzel, T.J., *Tetrahedron: Asymmetry*, **2010**, *21*, 2289-2294.

Association of Secondary Amines





In methanol



Dimethylbenzylamine


The *C*-methyl and *N*-methyl resonances (400 MHz) of (a) **3** (10 mM) with increasing concentrations of **1** (0, 5, and 10 mM), (b) the hydrochloride salt of **3** (10 mM) with increasing concentrations of **1** (0, 20, and 40 mM), and (c) **3** (10 mM) with increasing concentrations of L-tartaric acid (0, 5, and 10 mM)






^1H NMR spectrum (400 MHz) of the methine resonance of **8** (10 mM) in methanol- d_4 with **1** at (b) 5 mM, (c) 10 mM, (d) 15 mM, (e) 20 mM, (f) 30 mM, (g) 40 mM.

Lovely, A.E.; Wenzel, T.J., *Organic Letters*, 2006, 8, 2823-2826.




Pyrrolidines



Lovely, A.E.; Wenzel, T.J., *Tetrahedron Asymmetry*, **2006**, *17*, 2642-48











Piperidines and Piperazines



*Lovely, A.E.; Wenzel, T.J.,
Journal of Organic Chemistry,
2006, 71, 9178-82.*

Tertiary Amines





^1H – Discrimination usually small


^{13}C – Baseline discrimination

Cyclodextrins

- Cyclic oligosaccharides
- Glucose units
 - 6 – α
 - 7 – β
 - 8 – γ
- Water-soluble



Carboxymethylated Cyclodextrins Synthetic Schemes



Dignam, C.F.; Randall, L.A.; Blacken, R.D.; Cunningham, P.R.; Lester, S.-K.G.; Brown, M.J.; French, S.C.; Aniagyei, S.E.; Wenzel, T.J., *Tetrahedron Asymmetry*, **2006**, *17*, 1199-1208.



Wenzel, T. J.; Amonoo, E. P.; Shariff, S. S.; Aniagyei, S. E., *Tetrahedron: Asymmetry*, **2003**, *14*, 3099-3104.

Smith, K. J.; Wilcox, J. D.; Mirick, G. E.; Wacker, L. S.; Ryan, N. S.; Vensel, D. A.; Readling, R.; Domush, H. L.; Amonoo, E. P.; Shariff, S. S.; Wenzel, T. J., *Chirality*, **2003**, *15*, S150-S158.

Degree of CM Substitution



	2-position	6-position	Indiscriminate
α -CD	3	2	7
β -CD	4	1	9
γ -CD	4.5	2	8

α -CDCM





^1H NMR (400 MHz, D_2O) of (a) 10 mM chlorpheniramine with 10 mM (b) α -CD, (c) α -CDCM-Ind, (d) α -CDCM-2 and (e) α -CDCM-6.

β -CDCM



^1H NMR (400 MHz, D₂O) of (a) 10 mM chlorpheniramine with 10 mM (b) β -CD, (c) β -CDCM-Ind, (d) β -CDCM-2 and (e) β -CDCM-6. Impurities marked by "x"





–β-methylphenethyl-
amine HCl (10 mM)

–β-CM-CD-Ind (20 mM)

–Yb(III) – (2-8 mM)

K.A. Provencher, T.J. Wenzel, *Tetrahedron Asymmetry*, **2008**, *19*, 1797-1803
 Provencher, K.A.; Weber, M.A.; Randall, L.A.; Cunningham, P.R.; Dignam,
 C.F.; Wenzel, T.J., *Chirality*, **2010**, *22*, 336-346.

Cationic Cyclodextrins





DS = 0.7, 1.1, 1.5, 3.0

DS = 1.5 the best


C.D. Chisholm, T.J. Wenzel,
Tetrahedron Asymmetry, in press

Tyrosine with α -CD-GTAC



¹H NMR (400 MHz, D₂O) (a) 10 mM tyrosine (L>D), with (b) 10 mM native α -CD, (c) 20 mM native α -CD, (d) 10 mM α -CD GTAC (DS = 1.5), and (e) 20 mM α -CD GTAC (DS = 1.5).

Synthesis of Calix[4]resorcinarenes




Yanagihara, R.; Tominaga, M.; Aoyama, Y. *J. Org. Chem.*, **1994**, *59*, 6865-6867.


Dignam, C. F.; Zopf, J. J.; Richards, C. J.; Wenzel, T. J., *Journal of Organic Chemistry*, **2005**, *70*, 8071-8078.



Dignam, C. F.; Richards, C. J.; Zopf, J. J.; Wacker, L. S.; Wenzel, T. J., *Organic Letters*, **2005**, *7*, 1773-1776.

SCR-Pro







Water-soluble






10 mM substrate with 2, 4, 6, 8, and 10 mM SCR-Pro

Geometries of Complexes







Association Constants

	(<i>R</i>)- enantiomer	(<i>S</i>)- enantiomer
1-Phenylethylamine HCl	68	97
1-(1-naphthyl)ethylamine HCl	361	595
Propranolol HCl	258	482
Tryptophan methyl ester HCl	59	113
Sodium tryptophan	67	40



^1H NMR spectra (400 MHz, D_2O) of doxylamine (10 mM), (a) at 23°C , (b) with SCR-Pro (2 mM) at 23°C , (c) with SCR-Pro (40 mM) at 23°C and (d) with SCR-Pro (40 mM) at 50°C .

Other Calix[4]resorcinarene Derivatives




C.M. O'Farrell, J.M. Chudomel, J.M. Collins, D.F. Dignam, T.J. Wenzel, *Journal of Organic Chemistry*, **2008**, *73*, 2843-2851.

C.M. O'Farrell, T.J. Wenzel, *Tetrahedron Asymmetry*, **2008**, *19*, 1790-1796.


C.M. O'Farrell, K.A. Hagan, T.J. Wenzel, *Chirality*, **2009**, *21*, 911-921.

Hagan, K.A.; O'Farrell, C.M.; Wenzel, T.J., *European Journal of Organic Chemistry*, **2009**, 4825-4832.




2-*tert*-butylamino-1-phenylethanol (10 mM) with
6, 8, and 10 mM of **SCR-t4L**



Our Latest Calix[4]resorcinarene Systems



α -methyl-L-proline
(α MP)

L-pipecolic acid
(LPA)







(a) 2/3-(*S*), 1/3-(*R*)
(10 mM)
(b) t3L (30 mM)
(c) α MP (30 mM)


(c)

(b)


(a)



(a) 2/3-*L*, 1/3-*D*
(10 mM)
(b) c4L (10 mM)
(c) α MP (8 mM)



- (a) *N*-benzyl- α -methylbenzylamine
(10 mM) – CH₃
- (b) 20 mM t3L
- (c) 20 mM α MP
- (d) 15 mM LPA



- (a) 20 mM LPA
- (b) 2/3-D, 1/3-L (10 mM)

Special Thanks To:

- Catherine Dignam
- Kayla Zomlefer
- Jan Collins
- Alice Brechting
- Celeste Morin
- Rebecca Miles
- Daphney Frederique
- Melissa Roan
- Jeff Troughton
- Beth Pond
- Amanda Colby
- Matt Bogyo
- Estelle Lebeau
- Lauren Randall
- Renee Blacken
- Patrick Cunningham
- Shawna-Kaye Lester
- Monique Brown
- Susie French
- Stella Aniagyei
- Edwin Amonoo
- Sonia Shariff
- Ngoc Pham
- Cora Chisholm
- Kristin Smith
- James Wilcox
- Gudrun Mirick
- Lee Wacker
- Nicole Ryan
- David Vensel
- Ginger Readling
- Hilary Domush
- Jason Zopf
- Chris Richards
- Bailey Freeman
- David Sek
- Jolene Thurston
- Ann Lovely
- Courtney O'Farrell
- Katelyn Provencher
- Madeline Weber
- Brad Wenzel
- Katie Hagan
- Thao Nguyen
- Chloe Bourne
- Rebecca Clark
- Ryan Rollo

And To:


- National Science Foundation
- Research Corporation
- Camille and Henry Dreyfus Foundation
- Howard Hughes Medical Institute/Bates College
- Pfizer Pharmaceutical



THREE YEARS
OF $\Delta\Delta\delta\delta$



Pens and pencils are chiral (right-handed) because of the writing



Be careful what you put on a pencil!

